



## PICO Toner

### The hand held microcurrent device for home use

The updated PICO Toner has more available power to accommodate additional ports for additional applicators.

Now it has all-new accessory cables that are easier to connect and remove.

We did research and testing to come up with new advanced programming that offers updated signaling for better results!

The all-new PICO Toner has dual ports for adding accessories, up to 8 Alligator clips can be added or up to 8 Facial pads can be added.

#### Your PICO Toner Device Kit Includes

- x1 PICO Toner Device
- x1 Accessory Port Cable
- x1 Pair Silver Sculpting Gloves™ (one size fits most)
- x100 Poly Glove Inserts, (50 pairs)
- 4oz bottle of Ion Spray
- x1 Battery, AAA
- Wrist Strap for Hands Free Use
- Small Draw String Carrying Bag
- Online Training Videos and User's Quick Start Guide
- Customer Support & 1 Year Limited Warranty

01233 658 333

enquiries@neurotris.uk



## PICO Toner For Personal Facial & Body Toning

By popular demand, **Neurotris** introduces the world's most advanced personal face and body sculptor. The **Neurotris PICO Toner™** is the only hand held microcurrent device intended for home use featuring 2nd generation microcurrent "Constant Waveform Morphology" utilizing our patented Silver Sculpting Gloves™. If you are not ready for surgical procedures because of time, financial or medical limitations, we can offer a non-invasive approach to improve and return a youthful appearance. The experience is wonderfully relaxing and euphoric, causes no pain, swelling or downtime, and noticeable results begin immediately.

### The Most Versatile Personal Toner On The Market

The **PICO Toner™** is the only personal toner on the market today that allows for multiple output accessory attachments. The Pico Toner comes standard with the highest quality Silver Toning gloves available today. You have the option to add toning probes and or facial toning pads to personalize your treatment needs.

The toning gloves allows for more surface area to be treated as the probes can be used to treat fine areas. Use the Facial toning pads as a static treatment to focus on areas over the eye or nasolabial folds.

Not all microcurrent devices are considered the same. The advanced technology of the **Pico Toner** is dependent on the sophistication of engineering that takes place between the battery and the output accessories used (Gloves, Probes or Pads).

### Frequency Specific Preset Programs

The **PICO Toner™** For Personal Facial & Body Toning has two 30-minute programs that automatically alternate through THREE preset signatures to bring back the natural "glow" of youthful hues. No other service can come close to these proven results!

### Amplitude Control

The **PICO Toner™** For Personal Facial & Body Toning has user adjustable frequency amplitude control buttons to adjust the intensity during the session for ultimate personal comfort.

### Battery Life

The **PICO Toner™** For Personal Facial & Body Toning uses one AAA battery and has a long battery life. The smart feature built in automatically shuts off the unit after the 20 minute session.

### No More Gel!

The **PICO Toner™** For Personal Facial & Body Toning does not need expensive gels or creams to operate.



**01233 658 333**  
**enquiries@neurotris.uk**